



Sun•Flex - MacIntosh Flexible Tip Blades

English Profile

Laryngoscopes

With Conventional Illumination

SUN-FLEX™ (FLEXIBLE TIP)

- Satin finish
- Extra bright reflector lamp
- English channel to help visualize the epiglottis
- Surgical stainless steel construction
- Precise control to elevate the epiglottis
- Tip has non-pinching design

Sun•Med is the only manufacturer to provide both conventional and fiber optic flexible laryngoscope blades. Another Sun•Med product to assist you with difficult intubation. Blades made of surgical stainless steel.



REPLACEMENT LAMPS

5-0234-14

NEW! Parabolic Reflector Lamp
Provides 40% More Light!



CAT#	SIZE	DESCRIPTION	A	B
5-5082-03	3	Medium Adult	132mm	25mm
5-5082-04	4	Large Adult	155mm	26mm

CAT#	DESCRIPTION
5-5082-42	Sun-Flex MacIntosh Kit

KIT INCLUDES:

- Sun-Flex blades sizes 3 and 4
- 1 medium chrome-plated handle (stainless steel available)
- 1 extra reflector lamp
- Foam carrying case